

REFERENCE: ADDITIONAL CARD

REF NO:.....

REFERENCE CARD

REF NO: 112

TITLE: Methodology REF: Glasspoole EA, Ericson RL, Davidson CL
SUB TITLE: Sample preparation Demineralisation of enamel in relation to the
TOPIC: preparation of tooth surface fluoride release of materials. Am J Dent 2001;
14:8-12

"The middle 2mm of the facial enamel of each tooth was isolated between strips of attached vinyl tape for the purpose of effective demineralization and then for further remineralization. The rest of the tooth surface will be covered by nail varnish."